



# Routinely Used Surgical Guides

The expression 'routinely used' in the title does not suggest that the article below describes a simple, daily procedure. I think no single treatment ought to be labelled like this since neither a doctor/dentist nor the patient who is subject of it should ever have that feeling.


Here follows a report on a surgery accomplished with several guides which took place for the attendants on 9 July 2008 in the dental office of DentomiX-Ray at Moszkva tér in Budapest.



Dr Csigi Tamás and Dr Kecskés István had the object to demonstrate to colleagues how simply operations in implant surgery, compared to early stages, can be carried out with the guide planning software and milling machine by Schütz Dental and an E

Woo Tech 3D CT device.

With help of the software site, angulation, location and depth for the implants can be pre-determined. The freedom of the 3D imaging enables the user to calculate earlier unclear charac-



ters such as bone density, height of occlusion, etc.

The program was controlled by Schütz experts, the background for dental prosthetics was delivered by Royal Dental Studio while the implant system (DIO) by Dental Market Professional. This product is available in the domestic market since the beginning of the year.

With help of modern technology visitor dentist colleagues could follow the events on a projector in the waiting room.

Procedures were routinely carried out and thanks to the guide treatment time was greatly shortened which resulted in significantly reduced pains for the patients.





As I learnt the e software system completed with the 3D implant modelling is currently accessible in only one Hungarian dental centre. Let me draw the attention of the readers to the opportunity that the 3D implant software developed by Schütz is compatible with every implant system available in Hungary today. The manufacturer's program includes 70 pieces of implant software in total and can be upgraded with further systems of your choice.

I recommend to interested readers to inform themselves about the benefits of the systems. It is worth doing it!

*Laczkó Tamás*



## Éljen a 3 dimenziós képalkotás szabadságával!

### IMPLANTÁCIÓS MODELLEZÉS, HOGY NE ÉRJE MEGLEPETÉS...!

#### Szolgáltatásaink:

- hagyományos és digitális röntgenkészítés
- 3 dimenziós CT készítése, melyet CD-n adunk át segédprogrammal együtt
- 3 dimenziós implantációs sablon készítése (rendszerfüggetlenül)

**Várjuk Önt vagy betegét Buda szívében,  
a Moszkva tér 4. alatt!  
Nyitva tartás: hétfőtől péntekig 8.30-19.00**



**E-WOO**  
TECHNOLOGY  
E-WOO Technology Co., Ltd.

### DentomiX-RAY 3D fogászati röntgen diagnosztikai centrum

1024 Budapest, Moszkva tér 4. fszt. 1. Tel.: 316-5087, mobil: 06-20-946-3463

E-mail: csigi@t-online.hu, web: www.dentomix.hu